

E-Pacer

The Rifton E-Pacer is the largest member of the Pacer family. Its sit-to-stand lift functionality removes a major barrier to gait training for large or highly dependent clients. The E-Pacer can be a life-changing solution for clients who have grown too large or too dependent to be safely transferred into manual gait trainers. Its strong and stable frame accommodates users up to 6' 5" and 350 lb.



Effortless adjustments

The lift column adjusts from 32¾" to 59". The legs expand to fit around large wheelchairs.

Transfers made easy

The E-Pacer's electric lift makes the initial transfer safe and simple even with large clients. Here is one possible sequence for a no-lift transfer:



Transfer the client onto the pelvic support in a seated transfer using the high straps.



Use the pelvic support for a sit-to-stand transfer.




Begin gait training.

Gait Tracker

The E-Pacer scale connects with a phone or tablet via Bluetooth to display scale data on Rifton's Gait Tracker app.

The Gait Tracker shows a running average of the weight measured by the scale and calculates average weight-bearing for a gait training session.



Mount a phone on the patient handgrip with a bicycle phone mount to make the scale visible to the client.

Components



Scale

Weigh your clients during routine transfers, or measure their weight-bearing during gait training. When the scale is paired with a tablet or phone via Bluetooth, track your client's weight-bearing progress with Rifton's Gait Tracker app.

Note: The scale must be included in the original purchase. It cannot be added later.



Odometer

No more counting tiles on the floor! The odometer makes it easy to record how far your client walks in feet or meters.



Front handle

The front handle allows a caregiver to maneuver and guide the E-Pacer. It includes an up/down switch to adjust the height of the body support.



Switch pole

The switch pole lets a caregiver operate the electric lift without bending down (when not using the front handle). Turn it around to give the client control of the lift.



Arm prompts

Versatile arm prompts include height, rotation, angle, lateral proximity and forward/backward adjustments to accommodate various positioning requirements.

| Adjustment range | Padded surface area |
|--|---------------------------------|
| (elbow to fist): small 9" – 12½" large 10½" – 15" | small 4½" x 8" large 5" x 9" |





Arm platforms

Arm platforms are a simple arm support alternative for clients who need less positioning. They are made of closed-cell molded foam for padded comfort and easy cleaning.

Padded surface dimensions: 13½" x 5½"



Hip positioner

The hip positioner is designed to encourage good positioning for gait training. It enables natural pelvic movement and lower extremity weight-shifting, and provides good abduction. A pad is available for added comfort.

Dimensions: see dimensions chart



Pelvic support

Available in three sizes, the pelvic support is used for sit-to-stand transfers and gait training. It can also be used for fall prevention, making it suitable for early mobility in recovery. It may fit the client best when the straps are crossed as shown in the photo.

Dimensions: see dimensions chart



Thigh straps

Thigh straps are used for seated transfers. They are machine washable and are available in narrow and wide styles.

Width: narrow 5", wide 7"



Ankle prompts

Ankle prompts attach to the bottom bar of the E-Pacer frame and have padded straps that can be adjusted to prevent scissoring. Movable clips on the ankle prompt bar regulate stride length and placement.




Accessories tote


A tough and attractive tote bag big enough (14" x 16" x 9") to keep all your E-Pacer components together when not in use. Made of black ballistic nylon with a zippered top, here's a bag that will hold up to whatever you throw in it and still look great.




E-Pacer dimensions

| User dimensions (inches) | | K660 |
|---|--|---------------|
| Girth | | 22 min–55 max |
|  | Key user dimension: girth | |
| | Important: User's weight or height must not exceed: 350 lb or 77 inches | |
| Item dimensions (inches) | | K660 E-Pacer |
| Overall length | | 44¾ |
| Overall width | | 31½–51¾ |
| Overall height | | 32¾–59 |
| E-Pacer frame weight (lb) | | 72 |
| Turning diameter | | 50 |
| Floor to top of base leg for furniture clearance | | 8¾ |
| Min. user armpit height (when standing or ambulating) | | 31½ |
| Max. user height (when standing or ambulating) | | 77 |
| Max. working load (lb) | | 350 |

Pelvic support

| User dimensions (inches) Indicated by the red arrow below | | | |
|---|-----------|------------|-----|
| Key user dimension: half of thigh circumference measured at the groin. | | | |
| K306 small | K316 med. | K326 large | |
| 6–10 | 8–14 | 10–18 | |
|  | | | |
| Product dimensions (inches) | | | |
| K306 small | K316 med. | K326 large | |
| a | 8½ | 12 | 13½ |
| b | 16 | 20 | 25 |
| c | 9 | 12 | 14 |

Hip positioner

| Product dimensions (inches)–without pad | | |
|--|------------|----|
|  | | |
| K667 small | K668 large | |
| a Crotch length | 2½ | 3½ |
| b Width | 1¼ | 1¾ |
| c Length | 9 | 11 |

The E-Pacer comes in your choice of five colors.



Red



Blue



Black



Lime



Pink

2025 order form

Use dimension chart to select appropriate size.

Scan for sizing and pricing:
<https://www.rifton.com/e-pacer>



| | | K660 |
|--|--------|---|
| * Frame includes one battery, one battery charger and charger power cord. (International power cords available on request.) | | K661 <input type="checkbox"/> |
| * Color | CHOOSE | Red <input type="checkbox"/> |
| | | Blue <input type="checkbox"/> |
| | | Black <input type="checkbox"/> |
| | | Lime <input type="checkbox"/> |
| | | Pink <input type="checkbox"/> |
| Scale (must be chosen at the time of the original purchase. It cannot be added later.) | | K662 <input type="checkbox"/> |
| Odometer | | K663 <input type="checkbox"/> |
| Front handle | | K665 <input type="checkbox"/> |
| Switch pole | | K666 <input type="checkbox"/> |
| Arm supports (pair) | CHOOSE | Small arm prompts <input type="checkbox"/> |
| | | Large arm prompts <input type="checkbox"/> |
| | | Arm platforms <input type="checkbox"/> |
| | | Arm platforms with handgrips <input type="checkbox"/> |
| Hip positioner | CHOOSE | Small <input type="checkbox"/> |
| | | Large <input type="checkbox"/> |
| Hip positioner pad | CHOOSE | Small <input type="checkbox"/> |
| | | Large <input type="checkbox"/> |
| Pelvic support | CHOOSE | Small <input type="checkbox"/> |
| | | Medium <input type="checkbox"/> |
| | | Large <input type="checkbox"/> |
| Thigh straps | CHOOSE | Narrow <input type="checkbox"/> |
| | | Wide <input type="checkbox"/> |
| Ankle prompts | | K684 <input type="checkbox"/> |
| Additional battery | | K307 <input type="checkbox"/> |
| Additional thigh straps | CHOOSE | Narrow single <input type="checkbox"/> |
| | | Narrow pair <input type="checkbox"/> |
| | | Wide single <input type="checkbox"/> |
| | | Wide pair <input type="checkbox"/> |
| Rifton accessories tote | | K522 <input type="checkbox"/> |

* Required