

Choose your Anchor



To order or get a quote:

Go online: <https://www.rifton.com/anchors-sizing>

Email: sales@rifton.com

Call: 800.571.8198

Users with low tone receive support and develop strength for improved functional sitting. Users with high tone or movement disorders gain stability for increased functional use of the free hand and arm. Anchors promote inclusion and social development for the user.

Anchors are short-term positioning aids and are not intended to be used for primary standing support or transfers.



Hand Anchor: The hand Anchor attaches to any smooth surface to achieve stability and positioning. No more drilling through wheelchair trays!

K821

Overall length: 5½"

Length of handhold with bulb: 4¼"

Diameter of handhold: 1¼"



Grab Bar: The horizontal grab bar attaches to any smooth, horizontal surface to assist in functional activities. It can be used to stabilize both hands at once.

K824

Overall length: 20"

Overall length with cups in: 13½"

Diameter of handhold: 1¼"

Length of bar: 10"



Elbow Anchor: The elbow Anchor provides stability to the forearm while leaving the hands free for functional activities. It can also reduce abnormal posturing of the shoulder girdle by maintaining the elbow position.

Right K815

Left K816

Overall length: 11"

Overall width: 7¼"

Width of arm trough: 4"

Length of arm trough base: 9½"



Wrist Anchor: For a more secure handhold and to help maintain hand placement, the wrist Anchor includes a wrist strap and pad.

K822

Overall length: 8"

Length of handhold with bulb: 4"

Diameter of handhold: 1¼"



Horizontal Anchor: The versatile horizontal Anchor is ideal for clients who can grasp more easily with their forearm pronated. Attach to any smooth table or tray surface for stability and positioning.

K823

Overall length: 9½"

Length of handhold with bulb: 4½"

Diameter of handhold: 1¼"



Arm Anchor: The padded arm Anchor stabilizes the forearm. The handgrip extends and rotates to any angle to meet multiple positioning needs.

Small right K811

Small left K812

Inside of thumb to elbow: 7½" – 12"

Overall length with handhold: 9" – 13"

Width of arm trough: 4"

Large right K813

Large left K814

Inside of thumb to elbow: 10" – 15"

Overall length with handhold: 11" – 16"

Width of arm trough: 4½"