

# 2024 order form

Use dimension chart to select appropriate size.




**To order or get a quote:**  
**Go online:** <https://www.rifton.com/e-pacer-sizing>  
**Email:** [sales@rifton.com](mailto:sales@rifton.com)  
**Call:** 800.571.8198




		K660	
<b>* Frame</b> includes one battery, one battery charger and charger power cord. (International power cords available on request.)		K661	\$5,290 <input type="checkbox"/>
<b>* Color</b>		CHOOSE	Red <input type="checkbox"/>
			Blue <input type="checkbox"/>
			Gray <input type="checkbox"/>
			Lime <input type="checkbox"/>
			Pink <input type="checkbox"/>
<b>Scale</b> (must be chosen at the time of the original purchase. It cannot be added later.)		K662	\$1,260 <input type="checkbox"/>
<b>Odometer</b>		K663	\$210 <input type="checkbox"/>
<b>Front handle</b>		K665	\$400 <input type="checkbox"/>
<b>Switch pole</b>		K666	\$265 <input type="checkbox"/>
<b>Arm supports</b> (pair)	CHOOSE	Small arm prompts	K305 \$560 <input type="checkbox"/>
		Large arm prompts	K315 \$640 <input type="checkbox"/>
		Arm platforms	K397 \$420 <input type="checkbox"/>
		Arm platforms with handgrips	K396 \$475 <input type="checkbox"/>
<b>Hip positioner</b>	CHOOSE	Small	K667 \$350 <input type="checkbox"/>
		Large	K668 \$385 <input type="checkbox"/>
<b>Hip positioner pad</b>	CHOOSE	Small	K541 \$115 <input type="checkbox"/>
		Large	K542 \$155 <input type="checkbox"/>
<b>Pelvic support</b>	CHOOSE	Small	K306 \$286 <input type="checkbox"/>
		Medium	K316 \$320 <input type="checkbox"/>
		Large	K326 \$340 <input type="checkbox"/>
<b>Thigh straps</b>	CHOOSE	Narrow	K308 \$184 <input type="checkbox"/>
		Wide	K338 \$184 <input type="checkbox"/>
<b>Ankle prompts</b>		K684	\$230 <input type="checkbox"/>
<b>Additional battery</b>		K307	\$167 <input type="checkbox"/>
<b>Additional thigh straps</b>	CHOOSE	Narrow single	K318 \$92 <input type="checkbox"/>
		Narrow pair	K308 \$184 <input type="checkbox"/>
		Wide single	K328 \$92 <input type="checkbox"/>
		Wide pair	K338 \$184 <input type="checkbox"/>
<b>Rifton accessories tote</b>		K522	\$151 <input type="checkbox"/>

\* Required

# E-Pacer dimensions

User dimensions (inches)		K660
Girth		22 min–55 max
	<b>Key user dimension: girth</b>	
	<b>Important:</b> User's weight or height must not exceed: 350 lb or 77 inches	
Item dimensions (inches)		K660 E-Pacer
Overall length		44¾
Overall width		31½–51¾
Overall height		32¾–59
E-Pacer frame weight (lb)		72
Turning diameter		50
Floor to top of base leg for furniture clearance		8¾
Min. user armpit height (when standing or ambulating)		31½
Max. user height (when standing or ambulating)		77
Max. working load (lb)		350

## Pelvic support

User dimensions (inches)			
Indicated by the red arrow below			
<b>Key user dimension:</b> half of thigh circumference measured at the groin.			
K306 small	K316 med.	K326 large	
6–10	8–14	10–18	
			
Product dimensions (inches)			
K306 small	K316 med.	K326 large	
<b>a</b>	8½	12	13½
<b>b</b>	16	20	25
<b>c</b>	9	12	14

## Hip positioner

Product dimensions (inches)–without pad		
	K667 small	K668 large
<b>a</b> Crotch length	2½	3½
<b>b</b> Width	1¼	1¾
<b>c</b> Length	9	11

The E-Pacer comes in your choice of five colors.



Red



Blue



Gray



Lime



Pink