

2024 order form for large

Use dimension chart to select appropriate size.



To order or get a quote:
Go online: <https://www.rifton.com/activity-sizing>
Email: sales@rifton.com
Call: 800.571.8198



		Standard R860 large		Hi/lo R870 large							
* Seat and back includes seatbelt (Hi/lo includes footboard and push handles)	CHOOSE	Adjustable	R861	\$760	<input type="checkbox"/>	R871	\$1,450	<input type="checkbox"/>			
		Adjustable w/ back spring	R862	\$880	<input type="checkbox"/>	R872	\$1,570	<input type="checkbox"/>			
* Base	CHOOSE	Without spring	R873	\$465	<input type="checkbox"/>	R883	\$3,140	<input type="checkbox"/>			
		With spring	R874	\$560	<input type="checkbox"/>	N/A					
* Arm supports (pair)	CHOOSE	Armrests	R868	\$270	<input type="checkbox"/>	R868	\$270	<input type="checkbox"/>			
		Forearm prompts	R879	\$700	<input type="checkbox"/>	R879	\$700	<input type="checkbox"/>			
* Adjustable legs	CHOOSE	Short	R808	\$75	<input type="checkbox"/>	N/A					
		Long	R894	\$85	<input type="checkbox"/>	N/A					
		Short with casters	R895	\$330	<input type="checkbox"/>	N/A					
		Long with casters	R896	\$360	<input type="checkbox"/>	N/A					
* Pads	CHOOSE	Tan	R875	\$285	<input type="checkbox"/>	R875	\$285	<input type="checkbox"/>			
		Blue	R876	\$285	<input type="checkbox"/>	R876	\$285	<input type="checkbox"/>			
		Red	R877	\$285	<input type="checkbox"/>	R877	\$285	<input type="checkbox"/>			
		Pink	R977	\$285	<input type="checkbox"/>	R977	\$285	<input type="checkbox"/>			
		Purple	R978	\$285	<input type="checkbox"/>	R978	\$285	<input type="checkbox"/>			
		Green	R979	\$285	<input type="checkbox"/>	R979	\$285	<input type="checkbox"/>			
Headrest	CHOOSE	Contoured	R897	\$260	<input type="checkbox"/>	R897	\$260	<input type="checkbox"/>			
		Flat	R898	\$260	<input type="checkbox"/>	R898	\$260	<input type="checkbox"/>			
		Adjustable winged	R881	\$290	<input type="checkbox"/>	R881	\$290	<input type="checkbox"/>			
Laterals (pair)	CHOOSE	Small	R828	\$220	<input type="checkbox"/>	R828	\$220	<input type="checkbox"/>			
		Large	R838	\$235	<input type="checkbox"/>	R838	\$235	<input type="checkbox"/>			
		Large with chest strap	R839	\$285	<input type="checkbox"/>	R839	\$285	<input type="checkbox"/>			
Push handles (recommended for all chairs with wheels)			R878	\$166	<input type="checkbox"/>	Included					
Tray (requires armrests)			R669	\$320	<input type="checkbox"/>	R669	\$320	<input type="checkbox"/>			
Handhold			R649	\$31	<input type="checkbox"/>	R649	\$31	<input type="checkbox"/>			
Chest strap (wide)			R807	\$105	<input type="checkbox"/>	R807	\$105	<input type="checkbox"/>			
Butterfly harness	CHOOSE	Regular	R889	\$165	<input type="checkbox"/>	R889	\$165	<input type="checkbox"/>			
		Slim cut	R826	\$165	<input type="checkbox"/>	R826	\$165	<input type="checkbox"/>			
Thigh belt			R813	\$110	<input type="checkbox"/>	R813	\$110	<input type="checkbox"/>			
Pelvic harness			R846	\$210	<input type="checkbox"/>	R846	\$210	<input type="checkbox"/>			
Hip guides (pair)			R867	\$230	<input type="checkbox"/>	R867	\$230	<input type="checkbox"/>			
Abductor	CHOOSE	Small	R884	\$140	<input type="checkbox"/>	R884	\$140	<input type="checkbox"/>			
		Large	R891	\$160	<input type="checkbox"/>	R891	\$160	<input type="checkbox"/>			
Adductors (pair)			R892	\$220	<input type="checkbox"/>	R892	\$220	<input type="checkbox"/>			
Leg prompt			R859	\$260	<input type="checkbox"/>	R859	\$260	<input type="checkbox"/>			
Footboard (recommended for chairs with long legs)			R857	\$540	<input type="checkbox"/>	Included					
Ankle straps (pair)			R855	\$105	<input type="checkbox"/>	R855	\$105	<input type="checkbox"/>			
Sandals (pair)			R964	\$120	<input type="checkbox"/>	R964	\$120	<input type="checkbox"/>			
Wedges (pair, require sandals)			R699	\$85	<input type="checkbox"/>	R699	\$85	<input type="checkbox"/>			
Lumbar and seat support kit			R888	\$99	<input type="checkbox"/>	R888	\$99	<input type="checkbox"/>			
Backrest filler pad			R856	\$135	<input type="checkbox"/>	R856	\$135	<input type="checkbox"/>			
Whitmyer headrest adapter			R804	\$63	<input type="checkbox"/>	R804	\$63	<input type="checkbox"/>			
Directional Assist Wheel retrofit kit			N/A			R965	\$250	<input type="checkbox"/>	R966 (5 kits)	\$1400	<input type="checkbox"/>

* Required



Room to grow

Large chair at its smallest dimensions

- 1 Backrest height 19"
- 2 Distance between laterals 9½"
- 3 Seat depth 15"
- 4 Seat width 11"
- 5 Armrest height above seat 7"
- 6 Seat height above footboard 11¼"



Large chair at its largest dimensions

- 1 Backrest height 24"
- 2 Distance between laterals 14"
- 3 Seat depth 20"
- 4 Seat width 17"
- 5 Armrest height above seat 10½"
- 6 Seat height above footboard 21¼"



Large standard base
R860 Rifton Activity Chair



Large hi/lo base
R870 Rifton Activity Chair

User dimensions (inches)

Height 50–74



Key user dimension: height

The user's overall height is a general guide to help you select the appropriate chair. Choose the model that allows for growth.

Important: Make sure that seat width, depth and height are adequate for each individual user, and that the user's weight does not exceed the maximum working load.

Item dimensions (inches)

R860 large Standard base

Frame width	short legs: 26¾ long legs: 26¾ short legs w/ casters: 27¾ long legs w/ casters: 27¾
Seat height above floor	short legs: 16½–20½ long legs: 20–24 short legs w/ casters: 16½–20½ long legs w/ casters: 20–24
Seat angle (tilt-in-space)–degrees	15° forward, 15° back
Backrest angle–degrees	5° forward, 20° back
Footboard knee angle–degrees	45°–110°
Seat height above footboard	11¼–21¼
Seat width with hip guides without hip guides	11–14 17
Seat depth	15–20
Armrest height above seat	7–10½
Distance between laterals	9½–14
Backrest height	19–24
Headrest height above seat	19½–29½
Max. working load (lb)	250

User dimensions (inches)

Height 50–74



Key user dimension: height

The user's overall height is a general guide to help you select the appropriate chair. Choose the model that allows for growth.

Important: Make sure that seat width, depth and height are adequate for each individual user, and that the user's weight does not exceed the maximum working load.

Item dimensions (inches)

R870 large Hi/lo base

Frame width	29¾
Seat height above floor	16–24
Seat angle (tilt-in-space)–degrees	15° forward, 25° back
Backrest angle–degrees	5° forward, 20° back
Footboard knee angle–degrees	45°–110°
Seat height above footboard	11¼–21¼
Seat width with hip guides without hip guides	11–14 17
Seat depth	15–20
Armrest height above seat	7–10½
Distance between laterals	9½–14
Backrest height	19–24
Headrest height above seat	19½–29½
Max. working load (lb)	250