

2025 order form


Use dimension chart to select appropriate size.

Scan for sizing and pricing:
<https://www.rifton.com/mobilestander>



The **conversion kit** converts the K110 mini stander to the K130 small. It includes the deck and complete body support of the K130. Request from customer service.

Conversion kit K116 \$850

User dimensions (inches)	K110 mini	K130 small	K150 med.	K170 large
Height	25–32	30–40	37–52	48–65
 Key user dimension: height Select the appropriate Mobile Stander by the user's overall height. Choose the model that allows for growth.	Important: User's weight must not exceed the maximum working load.			
Item dimensions (inches)	K110 mini	K130 small	K150 med.	K170 large
Overall length and width	25x25	25x25½	32x29	39x33½
Overall height	22½–28	23½–30	31–39	37½–49
Width without large wheels	19½	19½	24½	28¾
Wheel diameter	20	20	27	35
Width of body support	7–10	7–10	9–12	10–14
Height of body support	16½–22	22–29	29–38	36–48
Item weight (lb)	28	30	40	65
Max. working load (lb)	50	50	120	175



* **Mobile Stander** frame includes body support with padding, seat pad, straps and fleece.

Specify color:
 pink, blue or lime

		K110 mini	K130 small	K150 medium	K170 large
CHOOSE	K101	\$2,270 <input type="checkbox"/>	K103 \$2,560 <input type="checkbox"/>	K105 \$3,100 <input type="checkbox"/>	K107 \$3,700 <input type="checkbox"/>
	Pink	<input type="checkbox"/>	Pink <input type="checkbox"/>	Pink <input type="checkbox"/>	Pink <input type="checkbox"/>
	Blue	<input type="checkbox"/>	Blue <input type="checkbox"/>	Blue <input type="checkbox"/>	Blue <input type="checkbox"/>
	Lime	<input type="checkbox"/>	Lime <input type="checkbox"/>	Lime <input type="checkbox"/>	Lime <input type="checkbox"/>
Wheels (pair) (include brakes)	K113	\$1,355 <input type="checkbox"/>	K113 \$1,355 <input type="checkbox"/>	K115 \$1,540 <input type="checkbox"/>	K117 \$1,855 <input type="checkbox"/>
Mounting bracket	K133	\$75 <input type="checkbox"/>	K133 \$75 <input type="checkbox"/>	K173 \$90 <input type="checkbox"/>	K173 \$90 <input type="checkbox"/>
Arm prompts (pair) (require mounting bracket)	K512	\$580 <input type="checkbox"/>	K512 \$580 <input type="checkbox"/>	K512 \$580 <input type="checkbox"/>	K512 \$580 <input type="checkbox"/>
Hand loops (require mounting bracket)	Single	K578 \$100 <input type="checkbox"/>	K578 \$100 <input type="checkbox"/>	K578 \$100 <input type="checkbox"/>	K578 \$100 <input type="checkbox"/>
	Pair	K579 \$200 <input type="checkbox"/>	K579 \$200 <input type="checkbox"/>	K579 \$200 <input type="checkbox"/>	K579 \$200 <input type="checkbox"/>
Communication tray (requires mounting bracket)	K518	\$435 <input type="checkbox"/>	K518 \$435 <input type="checkbox"/>	K518 \$435 <input type="checkbox"/>	K518 \$435 <input type="checkbox"/>
Additional seat pad	Small	N/A	N/A	K135 \$165 <input type="checkbox"/>	K135 \$165 <input type="checkbox"/>
	Large	N/A	N/A	K155 \$175 <input type="checkbox"/>	K155 \$175 <input type="checkbox"/>
Sandals with wedges (pair)	Small	K114 \$125 <input type="checkbox"/>	K114 \$125 <input type="checkbox"/>	N/A	N/A
	Medium	N/A	N/A	K154 \$135 <input type="checkbox"/>	N/A
	Large	N/A	N/A	N/A	K174 \$145 <input type="checkbox"/>

* Required