

2024 order form

Use dimension chart to select appropriate size.




To order or get a quote:
Go online: <https://www.rifton.com/prone-sizing>
Email: sales@rifton.com
Call: 800.571.8198



* **Prone Stander** includes a pair of trunk laterals, hip strap, padding and footboard.

Specify color: red or blue

User dimensions (inches)	E940 med.	E950 large
Height	44 – 59	57 – 72
 Key user dimension: height Select the appropriate Prone Stander by the user's overall height. Choose the model that allows for growth. Important: User's weight must not exceed the maximum working load.		
Item dimensions (inches)	E940 med.	E950 large
Board length	30	39
Board width (top / bottom)	11 / 16	14 / 19
Top of board to footboard	33 – 43½	44 – 57
Height when horizontal	29½	31½
Height when vertical	45½	61
Distance between laterals	8 – 16	12 – 22
Base width	26	30
Base length	39½	44
Tray inside width (widest point)	22½	29½
Tray inside length	17½	17½
Max. working load (lb)	150	200

		E940 medium		E950 large	
CHOOSE		E941	\$2,780 <input type="checkbox"/>	E951	\$3,400 <input type="checkbox"/>
		Red	<input type="checkbox"/>	Red	<input type="checkbox"/>
		Blue	<input type="checkbox"/>	Blue	<input type="checkbox"/>
Additional trunk laterals (pair)		E614	\$240 <input type="checkbox"/>	E614	\$240 <input type="checkbox"/>
Knee laterals (pair)		E942	\$224 <input type="checkbox"/>	E942	\$224 <input type="checkbox"/>
Additional hip strap		E946	\$88 <input type="checkbox"/>	E956	\$100 <input type="checkbox"/>
Tray		E947	\$675 <input type="checkbox"/>	E957	\$765 <input type="checkbox"/>
Round abductor		E933	\$114 <input type="checkbox"/>	E933	\$114 <input type="checkbox"/>
Collar (requires round abductor)		E934	\$94 <input type="checkbox"/>	E934	\$94 <input type="checkbox"/>
Round abductor with hip stabilizer		E655	\$330 <input type="checkbox"/>	E655	\$330 <input type="checkbox"/>
Adjustable abduction wedge		E945	\$270 <input type="checkbox"/>	E945	\$270 <input type="checkbox"/>
Adjustable abduction wedge with hip stabilizer		E605	\$345 <input type="checkbox"/>	E605	\$345 <input type="checkbox"/>
Sandals (pair)	CHOOSE	Small	E422 \$104 <input type="checkbox"/>	N/A	
		Medium	E433 \$114 <input type="checkbox"/>	E433	\$114 <input type="checkbox"/>
		Large	N/A	E432	\$114 <input type="checkbox"/>
Wedges (pair, require sandals)	CHOOSE	Small	E424 \$73 <input type="checkbox"/>	N/A	
		Medium	E435 \$78 <input type="checkbox"/>	E435	\$78 <input type="checkbox"/>
		Large	N/A	E434	\$83 <input type="checkbox"/>
Hand Anchor		K821	\$172 <input type="checkbox"/>	K821	\$172 <input type="checkbox"/>

* Required