

2024 order form

Use dimension chart to select appropriate size.




To order or get a quote:
Go online: <https://www.rifton.com/tram-sizing>
Email: sales@rifton.com
Call: 800.571.8198




			K310 Standard <small>(better for gait training)</small>		K320 Low-base <small>(fits under low beds)</small>			
* Frame includes one battery, one battery charger, accessories tote and charger power cord. (International power cords available on request.)	CHOOSE	Frame with scale	K301	\$6,910	<input type="checkbox"/>	K321	\$6,910	<input type="checkbox"/>
		Frame	K302	\$5,650	<input type="checkbox"/>	K322	\$5,650	<input type="checkbox"/>
* Thigh straps (pair)	CHOOSE	Narrow	K308	\$184	<input type="checkbox"/>	K308	\$184	<input type="checkbox"/>
		Wide	K338	\$184	<input type="checkbox"/>	K338	\$184	<input type="checkbox"/>
Switch pole			K666	\$265	<input type="checkbox"/>	K666	\$265	<input type="checkbox"/>
Arm supports (pair)	CHOOSE	Small arm prompts	K305	\$560	<input type="checkbox"/>	K305	\$560	<input type="checkbox"/>
		Large arm prompts	K315	\$640	<input type="checkbox"/>	K315	\$640	<input type="checkbox"/>
		Arm platforms	K397	\$420	<input type="checkbox"/>	K397	\$420	<input type="checkbox"/>
		Arm platforms w/ handgrips	K396	\$475	<input type="checkbox"/>	K396	\$475	<input type="checkbox"/>
Pelvic support	CHOOSE	Small	K306	\$286	<input type="checkbox"/>	K306	\$286	<input type="checkbox"/>
		Medium	K316	\$320	<input type="checkbox"/>	K316	\$320	<input type="checkbox"/>
		Large	K326	\$340	<input type="checkbox"/>	K326	\$340	<input type="checkbox"/>
Hip positioner	CHOOSE	Small	K667	\$350	<input type="checkbox"/>	K667	\$350	<input type="checkbox"/>
		Large	K668	\$385	<input type="checkbox"/>	K668	\$385	<input type="checkbox"/>
Hip positioner pad	CHOOSE	Small	K541	\$115	<input type="checkbox"/>	K541	\$115	<input type="checkbox"/>
		Large	K542	\$155	<input type="checkbox"/>	K542	\$155	<input type="checkbox"/>
Additional battery (recommended for institutional use)			K307	\$167	<input type="checkbox"/>	K307	\$167	<input type="checkbox"/>
Additional thigh straps	CHOOSE	Narrow single	K318	\$92	<input type="checkbox"/>	K318	\$92	<input type="checkbox"/>
		Narrow pair	K308	\$184	<input type="checkbox"/>	K308	\$184	<input type="checkbox"/>
		Wide single	K328	\$92	<input type="checkbox"/>	K328	\$92	<input type="checkbox"/>
		Wide pair	K338	\$184	<input type="checkbox"/>	K338	\$184	<input type="checkbox"/>
Swivel locks (pair)			K319	\$212	<input type="checkbox"/>	K324	\$212	<input type="checkbox"/>

* Required


TRAM dimensions

User dimensions (inches)		K310 and K320	
Girth		22 min–55 max	
 <p>Key user dimension: girth</p> <p>Important: User's weight or height must not exceed 350 lb or 76 inches</p>			
Item dimensions (inches)	K310 TRAM	K320 TRAM	
Overall length	45½	Same	
Overall width	27½–40	Same	
Overall height	43½–58	Same	
TRAM weight (lb)	70	Same	
Turning diameter	50	Same	
Floor to top of base leg for furniture clearance	6¾	4¼	
Min. user armpit height (when standing or ambulating)	30½	Same	
Max. user height (when standing or ambulating)	76	Same	
Max. working load (lb)	350	Same	

Pelvic support

User dimensions (inches)			
Indicated by the red arrow below			
Key user dimension: half of thigh circumference measured at the groin.			
	K306 small	K316 med.	K326 large
	6–10	8–14	10–18
			
Product dimensions (inches)			
	K306 small	K316 med.	K326 large
a	8½	12	13½
b	16	20	25
c	9	12	14

Hip positioner

Product dimensions (inches)–without pad		
		
	K667 small	K668 large
a Crotch length	2½	3½
b Width	1¼	1¾
c Length	9	11